

# CARELINK™



NOW AVAILABLE FROM TASK,  
IRELAND'S FIRST DIGITAL IP UNIT  
FOR HOMES WITH BROADBAND & WIFI

For further details on any product  
please contact us on the number  
below. Or complete and detach  
this panel to post in a sealed, stamped  
addressed envelope, and one of  
our friendly staff will give you a call.

**Call** 01 8435889  
**UK Customers Call** 0845 304 5535



## In an Emergency help is now at hand...

**How would you cope with an accident,  
sudden illness or a fall?**

The CareLink™ button can save you in an  
emergency, in the event of an accident or during a  
moment of concern. Indeed in any situation when you  
cannot reach the phone, simply press it to summon  
help, 24 hours a day, seven days a week.

This makes it ideal for those living or working alone,  
or for those with a disability and for anyone wanting  
peace of mind.

The CareLink™ Social Alarm comes with a slim  
pendant button. It can be worn around the neck, as a  
wristband, or clipped to clothing. The pendant  
can operate in your home and garden.

The latest CareLink™ alarm operates with or  
without a landline, by having options for switching  
between landline, Internet Protocol (IP) or mobile  
cellular communications (GSM).

# CARELINK™



## TASK S Community Care

TASK Community Care  
Unit 1, Block 3,  
City North Business Park  
Stamullen, Co.Meath, Ireland

**Call** 01 8435889  
**Email** [info@taskltd.ie](mailto:info@taskltd.ie)  
**Web** [www.taskltd.ie](http://www.taskltd.ie)

## TASK S Community Care

**UK Customers Contact:**  
TASK Community Care  
Catalyst Inc  
The Innovation Centre  
Queen's Road  
Belfast, BT3 9DT

**Call** 0845 304 5535  
**Email** [info@taskltd.com](mailto:info@taskltd.com)  
**Web** [www.taskltd.com](http://www.taskltd.com)

## Peace of mind & personal security at the touch of a button... with CareLink™

- Help at the touch of a button
- Easy to use
- Monitored 24 hours, 365 days a year
- Ideal for anyone of any age who feels at risk or vulnerable
- Affordable & reliable

## How your CareLink™ alarm works...

1. When you press the button, or when a detector is triggered, within seconds a call is placed to our CareLink Response Centre.
2. A two way voice channel is opened between the base unit in your home and our CareLink™ staff, twenty four hours a day, seven days a week.
3. CareLink staff will talk to you. If you tell them that something is wrong, or if they do not hear a response, they will immediately contact a family member or the emergency services, depending on the situation.



The CareLink™ Wristband can also be worn as a pendant or clipped to your belt. When pressed it transmits a wireless alert to the CareLink™ Base Unit which dials through to friendly operators at our 24 hour Monitoring Centre.



CareLink™ Wristband



CareLink™ Base Unit

Many additional telecare devices (peripherals) can be added wirelessly to your CareLink™ base unit, as and when you need them, giving your home additional protection. All these devices will send an alert through to our 24/7 monitoring centre when activated, turning any home into a 'smart home'. Our wide range of wireless telecare devices include those below:



## Any Questions?

TASK Community Care continues to make a major contribution to affordable community care for all, from our 'CareLink™ Two Way Speech Unit' to our many other telecare systems, all of which can be connected to our 24 hour monitoring centre.

Our telecare systems are used by individuals, local authorities, health professionals and in sheltered accommodation throughout Ireland & the UK, with acknowledged expertise and product reliability.

For further information on any of our products please contact us. Just place this slip in a stamped addressed envelope and send to your nearest TASK office.

Name

Address

Postcode

Telephone

Fax

Email

**CARELINK™**